

XRI®

X-Ray Inc. of Rhode Island

To place an order or inquire about pricing,
Please call: **(800) 556-6635 (401) 732-4300**

FAX: **(800) 543-1816**

or

E-mail: accessories@xrin.com

X-Ray Grid Encasements

Snap-On X-ray Grid Encasement (Grid Cap)

- Steel Construction with Lexan encased x-ray grid.
- Performs dual function of protecting X-ray grid with converting existing film cassette into a grid cassette.
- X-ray Cassette is held in place by pressure springs.
- Available in all sizes to fit all cassettes up 14 x 17".



Flat X-ray Encasement no Flange

Soyee X-ray grids may be installed in Flat Lexan encasements, which can be laid on a cassette or installed in the grid channels of a wall holder.

- Also available with flanged edges to fit wall holders which have a fixed 1/4" grid channel. .



XRAC X-ray Grid Encasement with Channels

- Soyee X-ray Grid is completely encased in White Lexan® with Polyethylene Channels.
- Allows you to easily slide the x-ray cassette into the channels while offering advanced protection to the Soyee X-ray grid.
- Available 8x10, 10x12, 24x30, 11x14 & 14x17.



XRAC EZ - SLIDE II

- Constructed of durable light weight aluminum with closed welded corners.
- Offers protection to your valuable Soyee X-ray grid while easily allowing the cassette to slide in and out.
- Resilient for Portable and Mobile units.
- Available 8x10, 10x12, 24x30, 11x14 and 14x17.



XRAC EZ - SLIDE III

- Constructed of durable light weight aluminum with closed welded corners.
- Fits into all standard portable x-ray cassette storage drawers.
- Available 14 x 17



XRI®

X-Ray Inc. of Rhode Island

To place an order or inquire about pricing,
Please call: **(800) 556-6635 (401) 732-4300**
FAX: **(800) 543-1816**

or

E-mail: accessories@xrin.com

Visit us at: <http://www.xriaccessories.com/>