



X-Ray Inc. of Rhode Island

To place an order or inquire about pricing,  
 Please call: **(401) 732-4300 (800) 556-6635**  
 FAX: **(401) 738-4708**

or

E-mail: [accessories@xrin.com](mailto:accessories@xrin.com)

**Call for final price**



Criteria	Deluxe Reclining	Deluxe Non-reclining
Seat Height	23" to 32"	23" to 32"
Mechanism	Hydraulic lift	Hydraulic lift
Arm Rests	Removable	Removable
Head Rest	Adjustable removable	N/A
Foot Rest	Removable	Removable
Chair Back	Reclines to 170 degrees	N/A
Chair Base	Five star, porcelain covered cast iron	Five star, porcelain covered cast iron
Castors	Locking, five	Locking, five
Foot Control	Raises/lowers patients, locks seat rotation	Raises/lowers patients, locks seat rotation
Color	Light gray	Light gray
Finish	High density powder coated	High density powder coated
Seat Cover	Fire retardant, stain resistant vinyl	Fire retardant, stain resistant vinyl
Shipping Weight	Approximately 95 lbs.	Approximately 95 lbs.
Weight Limit	300 lbs. maximum	300 lbs. maximum

**Mammography Chair Non-Reclining Model 23500**  
**Mammography Chair Reclining Model 23110**

- Designed for patient comfort and to satisfy the needs of the practitioner.
- Provided with a tubular footrest, polyurethane armrests and black vinyl upholstery.
- A hydraulic pump, controlled by a foot treadle, which is accessible from either side or back of the chair, provides a 6 1/2" adjustment from 21 1/2" to 28". The chair may be raised, lowered or locked anywhere within its 360° range of rotation.
- Easily removed flat headrest or a contoured headrest, which provides a full range of forward or backward movement, angular adjustment and a pivoting headrest cushion. Both headrests are height adjustable.
- Chair back may be reclined to 150° and to any intermediate position.
- Arm rests are removed by turning a knob that is under the seat. Once loosened, pull the armrest away from the chair. To put it back, slide it into place and tighten the knob. The arm rest can be removed prior to, or after, seating the patient.



**Model 23110 Shown**

**Standard in Black.**

Custom upholstered at additional cost, in choice of burgundy, deep brown, dove gray, dusty jade, gunmetal, navy blue, parchment, regimental blue, sea foam, space blue, tea rose, teal, violet and white. [Click here for Upholstery Chart.](#)

**Patient Weight Capacity is 250Lbs**

**Dimensions:**

Seat: 19" W x 19" D, Back: 18" W x 16" H, Overall Width: 24"  
 Depth of chair: (with footrest, reclined): 50"

Mammography Chairs - Removable Armrest with Five 2" Locking Dual Wheel Casters		
No. 23500	Chair with non-reclining backrest	Without headrest
No. 23110	Chair with reclining backrest	With flat headrest



Precision Skin Marking Systems

## Mammography Jackets



- Pink border identifies as separate from master jacket.
- Size: 11" x 14", including  $\frac{3}{4}$ " scored self expansion.
- Pocket size:  $8\frac{1}{4}$ " x 11".
- Illustrated back allows for Dx & patient communication.
- Offset double thumbcut for easy access to contents.
- Open side (opens along upper edge).

cat.#	qty./ctn	wt./ctn	1-9 ctns	10-24 ctns	25+ ctns
JJM-PPM	200	30#	\$159.00	\$154.00	\$149.00

## Mammography Marker Sets



### Screening ACR Mammo Set

(standard set shown)

ITEM	Cat.#	PRICE
• Standard Set	ACR-525	\$160.00
• Velcro™ Set	ACR-526	\$175.00
• Suction Cup Set	ACR-527	\$190.00
• Magnetic Set*	ACR-528	\$185.00



### Basic ACR Mammo Set

(suction cup set shown)

ITEM	Cat.#	PRICE
• Standard Set	ACR-510	\$295.00
• Velcro™ Set	ACR-511	\$300.00
• Suction Cup Set	ACR-512	\$300.00
• Magnetic Set*	ACR-513	\$300.00

#### ACR STANDARDIZED ABBREVIATIONS

Sets are designed respective of the ACR abbreviations.

<b>PREFIX</b> (7/16" letters)	R (Right) L (Left) M (Magnification)
<b>BASE</b> (1/4" letters)	MLO (Mediolateral Oblique) CC (Craniocaudal) ML (Mediolateral) LM (Lateromedial) XCC (Exaggerated Craniocaudal) CV (Cleavage) AT (Axillary Tail) TAN (Tangential) ID (Implant Displaced) SIO (Superolateral to Inferomedial Oblique) LMO (Lateromedial Oblique) FB (From Below)
<b>SUFFIX</b> (1/4" letters)	RL (Roll Lateral) RM (Roll Medial) RS (Roll Superior) RI (Roll Inferior)
<b>DEGREES</b> (1/4" letters)	(no standard)
<b>INITIALS</b> (1/4" letters)	(no standard)

\*Magnetic Sets include markers with magnetic backing that stick to a white metal wallboard which is attached with Velcro™ to the wall or equipment. Can be mounted horizontally or vertically.



**NEW COLOR!**  
**Mammography Pink!**


## Mammography Markers custom or replacement

Any marker can be created with tech's initials!

CPT-67 with initials



description	example	STANDARD		VELCRO™		SUCTION CUP		MAGNETIC*	
		cat.#	price	cat.#	price	cat.#	price	cat.#	price
single letter	R	CPT-01	\$13.00	CPT-01V	\$14.00	CPT-01S	\$15.00	CPT-01M	\$14.00
2-3 letters	R ML	CPT-23	\$14.00	CPT-23V	\$15.00	CPT-23S	\$16.00	CPT-23M	\$15.00
4-5 letters	L 50°	CPT-45	\$17.00	CPT-45V	\$18.00	CPT-45S	\$19.00	CPT-45M	\$18.00
6-7 letters	R CC MBT	CPT-67	\$20.00	CPT-67V	\$21.00	CPT-67S	\$22.00	CPT-67M	\$21.00
8+ letters	L XCC MBT	CPT-80	\$24.00	CPT-80V	\$25.00	CPT-80S	\$26.00	CPT-80M	\$25.00



- Same style as those on the ACR sets shown above.
- Aluminum-backed, high density letters.

- R & L are 7/16" - all other letters are 1/4".
- Comes with a hole for hanging (except magnetic version).

## Economy Nipple-Artifact Markers

Cat.# LNA-20  
\$50.00



- Encapsulated 2mm LEAD ball.
- Used to mark nipple, lesion, scar, artifact, etc.
- Neutral color, peel-off, self-adhesive backing.
- Non-sterile.
- 100 markers mounted on ten convenient strips.

- Maintains patient modesty and confidence.
- Comfortable, non-irritating materials.

**REDUCED PRICE!**

- Provides protection against cross-contamination.
- Eliminates laundry and replacement costs.



*Cloth-like comfort and performance provides security and dignity for your patients.*



### Mammo Exam Vest

Cat.#GMV-30  
\$79.00/case **NEW!**

- Vest style with open front for easy on/off.
- Snap closure front is quick and secure.
- Multi-layer, non-woven material for softness and strength.
- Size: 36" x 30".
- Color: BOLD BLUE.
- 50 vests per case.

### Modesty Smock

Cat.#GMS-27  
\$61.00/case **NEW!**

- Slip-over poncho style with open sides.
- Elastic top keeps smock comfortable and secure.
- Multi-layer, non-woven material for softness and strength.
- Size: one size fits most.
- Color: BOLD BLUE.
- 25 smocks per case.



### Mammo Exam Cape

Cat.#GMC-21  
\$39.00/case

- Open front/back.
- Soft 3-ply tissue.
- Size: 30" x 21".
- Color: WHITE.
- 100 capes per case.

### Modesty Drape

Cat.#GMD-26  
\$42.00/case **NEW!**

- Provides modesty for the upper body.
- Soft 2-ply tissue.
- Size: 24" x 40".
- Color: WHITE.
- 100 drapes per case.



### Mammo Exam Cape

Cat.#GMC-22  
\$46.00/case

- Open front/back.
- Soft 3-ply tissue.
- Size: 30" x 21".
- Color: PINK.
- 100 capes per case.





## M-Spray™

Cat.#PDM-509  
\$114.00 / case

- Designed to Clean and Disinfect **Mammography** equip.
- Contains NO alcohol, phenol or glutaraldehyde.
- Effective bactericide, fungicide, & virucide.
- Use after each patient.
- For compression plates, cassettes, etc.
- Kills MRSA, Avian Flu & HIV-1.
- Concentrated; just add water.
- 250ml bottles (1-dozen per case).



Mammographic Accreditation Phantom



Sony UPD-F750 Mammography

# WOMEN'S HEALTH IMAGING

*Integration of Mammography Priors*

**CADPRO<sup>®</sup>**  
Advantage



**XRI<sup>®</sup>**

X-Ray Inc. of Rhode Island

To place an order or inquire about pricing,  
Please call: **(401) 732-4300 (800) 556-6635**  
FAX: **(401) 738-4708**

or

E-mail: [accessories@xrin.com](mailto:accessories@xrin.com)

Visit us at <http://www.xriaccessories.com>